

D-GLUCONIC ACID

Product #: **GLA-F20 (20 Tests)**
GLA-F50 (50 Tests)

Enzymatic, UV-Method

INTENDED USE

D-Gluconic Acid FLEX-REAGENT™ is intended for measuring D-Gluconic acid, as a marker for Botrytis in infected grape juice and must.

KIT CONTENTS	20T	50T
Sample Blank (Buffer)	20 mL	50 mL
Chromogen Diluent	18 mL	45 mL
Chromogen	2 mL	5 mL
Enzyme: GK, 6PGDH	1 mL	2.5mL
D-Gluconic Acid Std, 2.0 g/L	2 mL	2 mL

SYSTEM REQUIREMENT

Spectrophotometer should be capable of reading 340 nm absorbance over a 0-2 A range with a 1 cm lightpath.

SAMPLES

Turbid sample should be centrifuged.

REAGENTS AND STORAGE

Kit contents are stable through the labeled expiration date when stored at 2-8 °C. Reagents are liquid and ready to use.

ASSAY PREPARATION

Prepare Working Reagent

Prepare sufficient Working Reagent (W-Rgt) as shown in the table, mix gently. W-Rgt is stable 4-weeks refrigerated, bring to room temperature before use.

Working Reagent	# of Cuvettes		
	10	30	50
Chromogen Diluent	9	27	45
Chromogen	1	3	5

TESTING PROCEDURE

- Two Zero-Rgt Cuvettes will be used to zero the spec. Label a Zero-Rgt Blank & Reaction Cuvette. Pipet DI Water into each of these cuvettes, as shown in table.

ADD	Volume/Cuvette			
	Zero-Rgt Cuvettes		Standard(s) & Sample(s)	
	Blank Cuvette	Reaction Cuvette	Blank Cuvette	Reaction Cuvette
DI Water	10µL	10µL	-	-
Std or Sample			10µL	10µL
Sample Blank	1000µL	-	1000µL	-
W-Rgt	-	950µL	-	950µL
Enzyme		50µL		50µL
Mix, wait 10-15 minutes Zero Spec & Read Blank Cuvettes - refer to Step #4 Zero Spec & Read Reaction Cuvettes - refer to Step #5				
Zero Spec & Read ABS	A = 0 "Zero Rgt Blank Cuvette"	A = 0 "Zero Reaction Cuvette"	A _{SAMPLE-BL}	A _{SAMPLE}

- For each Standard and each Sample, label a Blank & Reaction Cuvette. Pipet each Standard & Sample into both (BL & Rx) sets of cuvettes, as shown in the table.
- Dispense Sample Blank and W-Rgt to respective cuvettes (refer to Table) and incubate at 37C.
- Zero the Spec with **Zero-Rgt Blank Cuvette**; read Standard & Sample Blank cuvettes.
- Re-zero Spec with **Zero Reaction Cuvette**; read Standard and Sample Reaction cuvettes.

CALCULATIONS

Calculate Net ABS values by subtracting corresponding Blank ABS values from Reaction ABS values:

$$\text{Net } A_{\text{STD}} = A_{\text{STD}} - A_{\text{STD-BL}}$$

$$\text{Net } A_{\text{SAMPLE}} = A_{\text{SAMPLE}} - A_{\text{SAMPLE-BL}}$$

Calculate Concentration (based on 2.0 G/L Standard)

$$\text{D-Gluconic Acid G/L} = (2.0) \times (\text{d.f.}) \times \frac{\text{Net } A_{\text{SAMPLE}}}{\text{Net } A_{\text{STANDARD}}}$$

LINEARITY: For samples with initial results > 2 g/L, dilute and re-test the sample, multiply result by the dilution factor.

APPENDIX

TESTING PROCEDURE (Automation)

Contact Unitech Scientific for the ChemWell for Wine™ automated test procedures and technical support.

CHEMICAL PRINCIPLES

In presence of ATP and GluK (gluconate kinase), D-gluconic acid is phosphorylated to form D-Gluconate-6-Phosphate, which is subsequently decarboxylated by NADP in the presence of 6PGDH (6-P-gluconate-dehydrogenase). The reaction product NADPH is proportional to the concentration of D-Gluconic acid in samples, and is measured by the 340nm absorbance of the reaction.

MANUFACTURED BY: **UNITECH SCIENTIFIC**
 12026 Centralia Road Suite H, Hawaiian Gardens, CA 90716
 Tel: 562-924-5150 Fax: 562-809-3140
www.unitechscientific.com