

D-GLUCOSE

Hexokinase UV-Method

Product #:
GLU-F60 (30 Tests)
GLU-F150 (75 Tests)
GLU-F500 (250 Tests)

INTENDED USE

D-Glucose FLEX-REAGENT™ is intended for measuring D-glucose concentration in wine.

KIT CONTENTS

	<u>30T</u>	<u>75T</u>	<u>250T</u>
G/F Buffer	20 mL	50 mL	170 mL
HK/G6PDH Solution	0.4 mL	1 mL	3X1 mL
D-Glucose Std., 3 G/L	1 mL	5 mL	2X5 mL

SYSTEM REQUIREMENT

Spectrophotometer should be capable of reading 340 nm absorbance over a 0-2 A range with a 1 cm lightpath.

SAMPLES

If wine samples are visually clear, no sample pretreatment is needed. Filter or centrifuge turbid samples, e.g. juice, must or fermentation samples.

REAGENTS

Kit contents are ready to use are stable through the labeled expiration date when stored at 2-8 °C.

ASSAY PREPARATION

Working Reagent, manual:

(see Appendix for ChemWell for Wine™ **Working Reagent**)
 Prepare Working Reagent just prior to testing, based on the number of blank, standards and wine samples in your assay.

	<u>Per test</u>
G/F Buffer (#1)	0.67 mL
Deionized Water (=D.I.)	<u>1.33 mL</u>
Appx. Total Vol.	2 mL

Working reagent is stable for 2-days when refrigerated; allow it to reach room temperature prior to assay.

TESTING PROCEDURE

Pipet each solution (#1-4) into the cuvettes, as shown:

	<u>Blank</u>	<u>Standard</u>	<u>Sample(s)</u>
1. D.I. Water	10µL		
2. Standard/Sample		10µL	10µL
3. Working Reagent	2.0mL*	2.0mL	2.0mL

Mix cuvettes, incubate 3 minutes.

Zero spectrophotometer (340 nm) with Reagent Blank

Read A_{INITIAL} (Initial ABS)

*Note: 2.0mL = 2000µL

Mix HK/G6 Enzyme

4. HK-G6PDH Solution	10µL	10µL	10µL
-----------------------------	-------------	-------------	-------------

Mix cuvettes, incubate 20 minutes, Read A_{FINAL} (Final ABS).

CALCULATIONS

1. Our online "Flex Calculator™-GF" spreadsheet at <http://unitechscientific.com/calculators.htm> is available for download. G/L values will be calculated automatically.

2. Manual Calculation:

Calculate ΔA values and G/L as follows for each cuvette:

$$\Delta A = A_{\text{FINAL}} - A_{\text{INITIAL}}$$

Subtract the ΔA of the Reagent Blank from each sample and standard the ΔA:

$$\text{Net A} = \Delta A_{\text{SAMPLE}} - \Delta A_{\text{BLANK}}$$

Calculate GLU Concentration (based on **3.0 G/L** D-Glu Standard)

$$\text{D-Glu G/L} = (3.0) \times (\text{d.f.}) \times \frac{\text{Net A}_{\text{SAMPLE}}}{\text{Net A}_{\text{STANDARD}}}$$

The testing range is up to 8 G/L. If test result is over-range, dilute the sample with deionized (or distilled) water; re-assay & multiply this test result by the dilution factor.

QUALITY CONTROL

Test a sample (i.e. check wine or standard) with a known GF concentration in each assay to monitor assay performance. Reagents and technique are acceptable if results are within 15% of expected value. Factors that may affect the performance of this test include instrument function, temperature, glassware cleanliness, and pipetting inaccuracy.

APPENDIX

NOTES FOR ALTERNATE CALCULATIONS:

a. Extinction Coefficient (results based on factor; compare standard result with known value to verify recovery.)

$$\text{D-Glu (G/L)} = \text{Net A} \times 5.92$$

Factor is derived as follows:

$$\text{Glu (g/L)} = \frac{\text{Net A} \times \text{MW} \times \text{T.V.} \times \text{d.f.}}{(\epsilon)(P)(1000\text{mg/g})(\text{SV})}$$

Where:

MW =180.16G/mole

TV = total reaction volume (mL)

SV = sample volume (mL), See Procedure Step 2

ε (absorptivity of NADP) = 6.22 @334-340nm
 [or 3.4 @ 365nm]

P = 1 cm light path

d.f. = dilution factor (Refer to Note 3)

$$\text{D-Glu} = \frac{\text{Net A} \times 180.16 \times 2.02}{6.22 \times 1 \times 1000 \times 0.01} = 5.85 \times \text{Net A}$$

Sample volume inaccuracy will affect results with the extinction coefficient calculation method; use calibrated micropipettes.


b. Multi-point standard curve Sample concentrations are calculated from the best-fit standard curve. A 5-Level D-Glucose Standard kit is available from Unitech Scientific.

METHODOLOGY & CHEMICAL PRINCIPLES

Hexokinase (HK) catalyses the phosphorylation of D-glucose by adenosine-5'-triphosphate (ATP). Glucose 6 phosphate (G6P) and adenosine diphosphate (ADP) are products of these reactions.



In the presence of glucose-6-phosphate dehydrogenase (G6PDH), G6P is oxidized by nicotinamide-adenine dinucleotide phosphate (NADP); 6-phospho-gluconate (6PG) and NADPH are reaction products, as shown below:



The increase in NADPH concentration is measured at 340nm and is the basis for calculation of D-Glucose concentration in the sample.

SIGNIFICANCE OF MEASUREMENTS

Significance of Measurements: Reducing sugars are the predominant soluble components of soft fruits, with sucrose in low amounts. D-Glucose and D-fructose are the predominant reducing sugars in grape and other fruit juices. The ratio of glucose to fructose in mature grapes is "1", but ranges from 0.7-1.2 according to variety, maturity and fermentation conditions.

AUTOMATED TESTING

'ChemWell for Wine' Glu analysis is linear to 8 G/L and uses the following **Working Reagent**, which is stable for 8-hours refrigerated.

Placed the **Working Reagent** and HK/G6 Enzyme in CW reagent rack.

# of Tests	<u>25T</u>	<u>55T</u>	<u>90T</u>
G/F Buffer (#1)	3mL	5mL	8mL
Deionized Water	<u>5mL</u>	<u>8mL</u>	<u>13mL</u>
WRgt (Approx. Total)	8mL	13mL	21mL

(# of Tests accounts for Reagent Bottle dead volume)

Calculations:

'ChemWell for Wine' calculates results automatically from either one standard or a multi-point standard curve; dilutes and retests values above linear range.

TRADEMARKS:

"ChemWell for Wine", "Flex Calculator", "FLEX Reagent" are Trademarks of Unitech Scientific

MANUFACTURED BY: **UNITECH SCIENTIFIC**
12026 Centralia Road Suite H, Hawaiian Gardens, CA 90716
Tel: 562-924-5150 Fax: 562-809-3140
www.unitechscientific.com