



Y S I Environmental

YSI 5905-W or 5010-W Wine Bottle Oxygen Probe

Pure
Data for a
Healthy
Planet.™

To order or for more
information, contact
YSI Environmental.

800 897-4151

YSI Environmental
937 767 7241
Fax 937 767 9353
environmental@ysi.com

Endeco/YSI
508 748 0366
Fax 508 748 2543
environmental@ysi.com

YSI Hydrodata
European Support Centre
44 1462 673581
Fax 44 1462 673582
sales@hydrodata.co.uk

YSI (Hong Kong) Limited
852 2891 8154
Fax 852 2834 0034
hongkong@ysi.com

YSI/Nanotech (Japan)
81 44 222 0009
Fax 81 44 222 1102
nanotech@ysi.com

YSI (Qingdao) Limited
86 532 389 6648
Fax 86 532 389 6647
china@ysi.com

ISO 9001
ISO 14001

Who's Minding the Planet? and Pure
Data for a Healthy Planet are registered
trademarks of YSI Incorporated.

Printed in USA

YSI incorporated
Who's Minding
the Planet?®



- Uses screw-on cap membrane
- Fits any YSI laboratory DO meter or lab/field instruments (may need adapter)
- Probe length = 6 inches
- Probe body tapered to fit into the neck of a wine bottle
- Self-stirring probe – No stir bar needed!

YSI 5905-W (*shown*) wine bottle probe can be used with all grey-box style instruments and the YSI 5000/5100 lab-only instruments (5011 adapter required). YSI 5010-W is also available with a smaller connector that plugs directly to the YSI 5000/5100 lab-only instruments.

YSI 5905-W and YSI 5010-W Probe Specifications

Accuracy	mg/L	±0.1 mg/L
	% air	±1% of reading
	temperature	±0.2° C
Range	mg/L	0-20
	% air	0-200%
	temperature	0 to +45° C
Resolution	mg/L	0.01
	% air	0.1%
	temperature	0.1° C
Response	Typical response at 20°C is 90% of the change within 30 seconds.	

www.YSI.com